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## English

**INSTRUCTIONS FOR USE AND PRESERVATION**

Dear customer,

Thank you for your trust in this Orliman product. You have acquired a quality product with a high medicinal value. Please, read the instructions carefully. If you have any doubt, get in contact with your doctor, your orthopaedist specialist or our customer service department. Orliman appreciates your choice and wishes you a speedy recovery.

**REGULATION**

These products satisfy the European Medical Device Directive 93/42/CEE (RD. 1591/2009). A Risk Analysis has been carried out (UNE EN ISO 14971) minimising all the existing risks. The trials have been carried out in accordance with the European regulation UNE-EN ISO 22523 of Prosthesis and Orthosis.

**FITTING INSTRUCTIONS**

**THE PRODUCT SHOULD BE ADJUSTED BY YOUR ORTHOPAEDIC TECHNICIAN OR IN ACCORDANCE WITH DOCTOR'S ORDERS.**

In order to obtain a greater degree of therapeutic effectiveness in the different pathologies and to prolong the life and utility of the product, it is fundamental that the correct size adequate for each patient or user is chosen (enclosed in the packet is a diagram with all the sizes and the equivalent in centimetres). An excess compression can produce intolerance, therefore we recommend that the compression be regulated to an optimum degree.

The bracing antiequine Boxia, has two modes of adaptation.

A-Go directly to the shoe (external use)

B-band half foot (internal use).

It is advisable to put an interface in the foot (eg, a sock).

**Fixing the shoe-Ref.: AB01/AB02**

Attach the strap above the ankle suprameleolar matching gel pad with Achilles tendon. Open the shoelaces of the shoe and put the plate hook attachment between the tongue of the shoe and shoelaces, trying to stay focused and with the hook in the direction of traction. Then enter the stretch drive by the guide at the front of the ankle. Set the hook and proceed to pull the ends with velcro fixing it after obtaining the desired traction. Type of footwear laces.

**Fixing band-Ref.: AB12/AB22**

Attach the strap above the ankle suprameleolar matching gel pad with Achilles tendon. Set banding (left or right) in the middle of the foot velcro closures on the instep area, ensuring that the locking hook is in the right direction. Then enter the stretch drive by the guide at the front of the ankle. Set the hook and proceed to pull the ends with velcro fixing it after obtaining the desired traction.

**Fitting the Ref.: AB14/AB24**

1-Size selection.

2-The calf support must correspond with the size of the Boxia AB01/AB02, previously checked and adapted to fit the patient.

3-The calf support is fastened with Velcro to the inner and front surfaces of the Boxia supramalleolar support and needs to be aligned with the protective pad located in the instep area of the support.

4-The Y-shaped calf support consists of three straps. The short and the long narrow strap are wrapped and fastened around the top of the calf. The remaining strap passes down the front of the tibia and is fastened to the inner surface of the Boxia supramalleolar support, in alignment with the protective pad in the instep area.

**PRECAUTIONS**

The constructive material is inflammable.

Do not expose these products to situations where they could ignite. If this were the case remove them quickly and use the necessary measures to put it out. The material used is hypoallergenic, however we cannot guarantee 100% that in some determined cases it will not produce allergies, if this occurs, remove the product and see the doctor that prescribed it.

In case of minor irritations due to sweat, we recommend the use of a cotton sock, to separate the skin from contact with the material. For small inconvenience caused by sweat, we recommend using an interface to separate the skin from contact with the tissue. If rash, irritation or swelling to remove the product and consult a doctor or prosthetist. Contraindicated in open scars with swelling, redness and accumulation heat.

The product contains natural rubber latex and may elicit allergic responses in individuals who are sensitised to latex.

**RECOMMENDATION- WARNINGS**

These products are to be used only by those which have been medically prescribed to do so. Persons unconnected to such prescriptions should not use them. When the product is not being used, store in its original packaging. In order to dispose of the packaging and the product, strictly comply with the legal norms of your community.

The use of the orthosis is conditioned to the recommendations of the prescribing doctor, therefore it must not be used for any other means than those prescribed.

For the validity of the guarantee, the issuing establishment should comply with this sheet of instructions.

Includes metal components so extreme caution to other uses, and in case of Magnetic Resonance.

**MANUFACTURE - CHARACTERISTICS**

The material features used in production have been tested and approved and meet all the European quality requirements.

All the products are produced with high quality materials and offer an unbeatable comfort and quality of use. All products offer restraint, stability and compression for the optimum treatment of the pathologies for which they have been designed.

**RECOMMENDATIONS FOR PRESERVATION AND CLEANING**

Stick the velcros to each other (if the orthosis has any), wash periodically by hand with a damp cloth and a neutral soap.

To dry the orthosis, use a dry towel to absorb most of the dampness and leave to dry at room temperature. Do not hang up

or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sun light etc. During use or during washing do not use alcohols, ointments or dissolvent liquids. If the orthosis is not properly dried any detergent residues could irritate the skin and deteriorate the product.

**GUARANTEES**

ORLIMAN, S.L.U. Guarantees all its products, as long as these have not been manipulated or altered from their original state. Does not guarantee those products in which as a result of misuse, deficiencies or breakage of any kind, its original state has been altered. If you observe any deficiency or anomaly, communicate this immediately to the establishment from which it was obtained in order for it to be changed.

ORLIMAN, S.L.U. Appreciates your choice and wishes you a speedy recovery.

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