# 1H08A

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ISO 9001

# 1H08A

## Articulated foot with shock absorber and heel adjustment

### INSTRUCTIONS FOR USE AND MAINTENANCE

Dear Customer,  
Please read the instructions carefully. Make sure that the user understands the instructions for use, paying particular attention to the safety and maintenance warnings.

### APPLICATION


The 1H08A articulated foot has been designed for patients with unilateral or bilateral femoral amputation at different levels or tibial proximal, transtibial or distal amputation. It can also be used in prostheses after knee or hip disarticulation. The prosthesis is made for the right and left foot. The foot has two fixed heel height positions.

### TECHNICAL CHARACTERISTICS

Code	1H08A									
Foot dimensions	22	23	24	25	26	27	28	29	30	
Foot length, mm	220	230	240	250	260	270	280	290	300	
Articulated foot weight, g	527 ±10	557 ±10	587 ±10	617 ±10	667 ±10	697 ±10	717 ±10	757 ±10	817 ±10	
Rigidity of the fitted shock absorber	soft (S)					normal (N)				
For patients weighing	up to 100 kg									
Foot heel height (without shoes) on the supporting surface	Position 1 - 25 mm									
	Position 2 - 10 mm									
Temperature range for use	From -40 °C to +40 °C									
Angular adjustment of the whole prosthesis on the frontal and sagittal planes	not less than ±7.5°									
Warranty period	2 years									

### DESCRIPTION OF THE STRUCTURE

The 1H08A articulated foot consists (see Fig. 1) of support element 1, socket 2 with safety 12, limiter 9 and cosmetic cover 3 made from microcellular polyurethane, housing 7, joint 4, washer 5, nut 6, nut 11 and disc 10.

- |                  |                  |           |
|------------------|------------------|-----------|
| 1 Main body      | 5 Washer         | 9 Limiter |
| 2 Socket         | 6 Nut            | 10 Disc   |
| 3 Cosmetic cover | 7 Housing        | 11 Nut    |
| 4 Joint          | 8 Shock absorber | 12 Safety |

The kit for attaching the 2H04A/1 articulated ankle joint includes a set of rubber shock absorbers 8 with different degrees of rigidity (soft, normal and rigid), marked respectively with the letters S, N and R.

The articulated ankle joint is attached through the hole in the prosthetic foot and fixed into place using nut 6 and safety washer 5.

Prosthetic feet with sizes 22, 23 and 24 come with a soft (S) shock absorber 8; all other feet come with normal (N) rigidity shock absorbers 8. Should the customer request the replacement of shock absorber 8 with one of different rigidity (S, N or R), it will be supplied additionally. Ref 304242/2S (soft); Ref 304242/2H (hard).

The articulated foot is compatible with Orliman's modular lower-limb prosthetic systems.

### OPERATION

During the stance phase (forward impulse) of the gait cycle and under the weight of the body, shock absorber 8 is compressed and the forefoot drops onto the supporting surface. This ensures a more stable gait with the prosthesis. During the gradual transition towards the end of the stance phase (backward impulse), the deformation of the forefoot (tip of the foot) takes place.

To change the heel height, shock absorber 8 should be compressed by removing the load of limiter 9 and, through the foot's side holes, move limiter 9 to the stop in fixed position 1 or 2.

### ADJUSTMENT

Adjustment of the articulated foot is carried out by putting limiter 9 into position 1. Changing the rigidity of the heel is carried out by replacing shock absorbers 8 with others that have a different rigidity. Depending on the patient's weight, it is recommended to use the shock absorbers in the following way:

- up to 70 kg - soft (S) or normal (N)
- up to 100 kg - normal (N) or rigid (R)

To replace shock absorber 8, unscrew screw 11, replace the shock absorber, insert disc 10 and tighten nut 11 within the limits of dimension B.

It is also possible to adjust the rigidity of the heel. The degree of compression of shock absorber 8 can be adjusted by screwing or unscrewing screw 11 within the limits of dimension B.

### MAINTENANCE

During use or cleaning, do not use talcum powder. Do not put the product into contact with corrosive or abrasive elements such as fresh or salt water, acids, dust, sand or other substances.

At least once a year, have the product checked by an orthopaedic technician for proper operation and wear. If you detect any unusual changes in the product's functioning or hear it making noises, stop using it immediately and notify an orthopaedic technician.

### RECOMMENDATIONS - WARNINGS

These products should not be used for purposes other than those described in these instructions. They should only be used by one person and under the conditions prescribed.

Before use, the product must be adapted by an orthopaedic technician.

When disposing of the product and packaging, strictly comply with the refuse disposal regulations of your community.

In knee joints, there is risk of entrapment, use caution.

As the product features metal components, use caution in the case of magnetic resonance.

### WARRANTY

These products comply with European Medical Devices Directive 93/42/EEC (Spanish Royal Decree 1591/2009). They have undergone tests in accordance with international standard EN ISO 22523 for prostheses and orthoses and EN ISO 10328 for lower-limb prostheses.

ORLIMAN, S.L.U. guarantees all of its products, provided that they have not been tampered with or altered from their original form. The company does not guarantee products whose characteristics have been altered through any kind of misuse, breakage or defect. If you notice any defect or abnormality, notify the supplier immediately in order for the product to be replaced.

ORLIMAN S. L. U.  
C/ Ausias March, 3 - Pol. Ind. La Pobla-L'Eliaua - Apdo. de correos 49  
C.P.: 46185 - La Pobla de Vallbona Valencia - España (Spain)  
Tel.: +34 96 272 57 04 - Fax: +34 96 275 87 00 - Exports Tel.: +34 96 274 23 33  
E-mail: orto@orliman.com - Exports email: export@orliman.com - www.orliman.com

### WARRANTY

2 years from purchase

Reference	Batch no.	Date of purchase
<input type="text"/>	<input type="text"/>	<input type="text"/>
Signature and stamp of the orthopaedic establishment		For the warranty to be valid, it is necessary to complete these details.
<input type="text"/>		<input type="text"/>