4H16C - 4H16T

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4H16C - 4H16T

› Shift adapter

INSTRUCTIONS FOR USE AND MAINTENANCE

Dear Customer,
Please read the instructions carefully. Make sure that the user understands the instructions for use, paying particular attention to the safety and maintenance warnings

APPLICATION

The shift adapter is used in lower-limb prostheses and serves to connect and adjust the modules of the prosthesis. The shift adapter ensures linear movement and fastening of housing 1 with respect to housing 2 and enables connection and angular adjustment of the complementary modules.

TECHNICAL CHARACTERISTICS


Code	4H16C	4H16T
Material	Stainless steel	Titanium
Item weight, g	340	213
Patient weight, kg, no higher than	125	150
Length, mm	60	
Temperature range for use	De -40 a +40°C	
Warranty period	3 years	
Axis movement range of housing 1 with respect to housing 2 which supports the prosthesis (in the replacement of screw 5), mm	De -11° C a +11°C	
Angular adjustment of the whole prosthesis on the frontal and sagittal planes	De -7,5 a +7,5°	

DESCRIPTION OF THE STRUCTURE

The shift adapter (Fig. 1) consists of housing 1 with a bolt and housing 2 with a trapezoidal groove. Both are connected together by a male-female system secured with screws B. Housing 2 is movable in two directions with respect to housing 1 (see Fig. 2). Screw 5 is placed into the threaded hole of housing 1 in the direction of movement of housing 2.

Both housings have spherical support surfaces A and two pairs of screws 3 and 4. These connections enable angular adjustment of the prosthesis with other modular systems on the frontal and sagittal planes of -7.5 to +7.5.

The shift adapter is compatible with Orliman's modular lower-limb prosthetic systems.

- 1 Housing
- 2 Housing
- 3 M8 screw, L=12 mm
- 4 M8 screw, L=12 mm
- 5 M4 screw


FITTING THE ITEM TO THE PROSTHESIS

ATTENTION. When fitting the shift adapter to the prosthesis, fasten screws B to the outer side of the body.

After selecting the direction of movement of housing 2 with respect to housing 1, unscrew screw 5 and place it into the threaded hole of housing 1 in the direction of movement of housing 2. Apply a layer of Loctite 243 (or 241) threadlocker to the thread of screw 5 and tighten with a torque of 3 Nm.

Loosen screws B in the male-female connection to prevent locking. Unlock by gently tapping housing 1 with a small hammer from the side opposite screws B.

During connection through the spherical surface, the screws should be fixed according to Figure 3 as follows: the 12 mm screws should be placed where the space between the ring and the pyramid is smaller and the 16 mm screws where the space is greater, ensuring that the screws protrude from the outer surface of the double adapter.

a), b) - correct connection of the module

c) - incorrect connection of the module

ATTENTION. To avoid annoying friction noises, there should not be any contact between housings 1 and 2 of the adapter and the pyramid of the complementary modules.

During final assembly of the prosthesis, threadlocker to the threads of screws 3 and 4 and tighten with a torque of 15 Nm.

RECOMMENDATIONS - WARNINGS

These products should not be used for purposes other than those described in these instructions. They should only be used by one person and under the conditions prescribed.

Before use, the product must be adapted by an orthopaedic technician.

When disposing of the product and packaging, strictly comply with the refuse disposal regulations of your community.

In knee joints, there is risk of entrapment, use caution.

As the product features metal components, use caution in the case of magnetic resonance.

WARRANTY

These products comply with European Medical Devices Directive 93/42/EEC (Spanish Royal Decree 1591/2009). They have undergone tests in accordance with international standard EN ISO 22523 for prostheses and orthoses and EN ISO 10328 for lower-limb prostheses.

ORLIMAN, S.L.U. guarantees all of its products, provided that they have not been tampered with or altered from their original form. The company does not guarantee products whose characteristics have been altered through any kind of misuse, breakage or defect. If you notice any defect or abnormality, notify the supplier immediately in order for the product to be replaced.

ORLIMAN S. L. U.
C/ Ausias March, 3 - Pol. Ind. La Pobla-L'Eliaana - Apdo. de correos 49
C.P.: 46185 - La Pobla de Vallbona Valencia - España (Spain)
Tel.: +34 96 272 57 04 - Fax: +34 96 275 87 00 - Tel. Exportación: +34 96 274 23 33
E-mail: orto@orliman.com · Export mail: export@orliman.com · www.orliman.com

WARRANTY

3 years from purchase

Reference	Batch no.	Date of purchase
<input type="text"/>	<input type="text"/>	<input type="text"/>

Signature and stamp of the orthopaedic establishment

For the warranty to be valid, it is necessary to complete these details.

LB-I-4H16

EDITION: 05/2015

ISO 13485

ISO 9001