6H09A

› Rodilla policéntrica de 4 ejes

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6H09A

› Polycentric knee joint

INSTRUCTIONS FOR USE AND MAINTENANCE

Dear Customer,
Please read the instructions carefully. Make sure that the user understands the instructions for use, paying particular attention to the safety and maintenance warnings.

APPLICATION

The polycentric knee joint has been designed for femoral prostheses in cases of unilateral or bilateral amputation and is recommended for disabled people with low to medium activity levels.

TECHNICAL CHARACTERISTICS

Code	6H09A
Material of the main parts	Aluminium
Item weight, g	460
For patients weighing	up to 100 kg
Maximum angle of flexion	135°
Temperature range for use	from -40°C to +40°C
Estimated useful life	3 years
Enables angular adjustment on the frontal and sagittal planes (of the whole prosthesis)	
▶ for the upper support	in a range of -7.5° to +7.5°
▶ for the lower support	in a range of -5° to +5°
Warranty period	2 years


DESCRIPTION OF THE STRUCTURE

The polycentric knee joint (Fig. 1) consists of upper and lower supports 1 and 18 linked by connecting rod 2, a pair of connecting rods 17 and two pairs of axes 14 and 16. Upper support 1 and lower support 18 feature needle roller bearings 20.

Upper support 1 features a stop which includes regulating screw 12, installation screw 19 and stops 13. During the stance phase of the gait cycle, the stop is supported by shock absorber 11 of connecting rod 2. The stop serves to adjust patellar flexion stability.

Lower support 18 features the tibial hinging device: spring 4, one of whose ends is supported by regulating screw 5 screwed into bushing 3 and the other end, through ball 7 and rod 6, is supported by connecting rod 2. Regulating screw 5 enables adjustment of the compression force of spring 4.

The polycentric knee joint is compatible with Orliman's modular lower-limb prosthetic systems.

1 Upper support	8 Screw	2 unit	15 Screw	4 unit
2 Connecting rod	9 Pyramid		16 Axis	2 units
3 Bushing	10 Pyramid		17 Side connecting rods	2 units
4 Spring	11 Shock absorber		18 Lower support	
5 Regulating screw	12 Regulating screw		19 Installation screw	
6 Pin	13 Stops		20 Needle roller bearing	8 units
7 Ball	14 Axis	2 units		

OPERATION

During the stance phase of the gait cycle, in which the prosthesis is supported by the whole foot, shock absorber 11 makes contact with the gasket and the immediate centre of rotation is located above the knee's reference axis and behind the load application line, which provides a stable position for the disabled person on the sagittal plane.


At the beginning of the swing phase, knee flexion occurs and spring 4 of the tibial hinging device is compressed; the stored energy in spring 4 promotes the carrying forward of the tibia when the polycentric knee joint is straightened.

The kinematics of this module ensures the functional shortening of the prosthesis by up to 15 mm during the swing phase on the supporting surface, which enables the construction of a prosthesis with equal length to the healthy limb.

FITTING THE ITEM TO THE PROSTHESIS

The assembly diagram for the prosthesis using the polycentric knee joint is shown in Figure 2. The recommended height for heel H is shown in the prosthetic foot's accompanying documentation.

During assembly of the prosthesis and to ensure patellar flexion stability, with the possibility of subsequent adjustment, it is necessary to place the knee so that both lower axes are located on horizontal line A (see Fig. 2).

ADJUSTMENT

To adjust patellar flexion stability, installation screw 19 must be removed. When regulating screw 12 is turned clockwise using a hex key, patellar flexion stability decreases, and when turned anti-clockwise, it increases. After adjusting patellar flexion stability, it is necessary to lock regulating screw 12 with installation screw 19 by applying Loctite 243 threadlocker to the thread. Torque for installation screw 19 – (6±2) Nm.

Regulating the carrying forward speed of the tibia is carried out by changing the length of spring 4 – it is necessary to tighten or loosen screw 5 to a depth of 7 mm with respect to the lower end of bushing 3.

MAINTENANCE

During the period of estimated useful life, the knee joint does not require prophylactic maintenance.

In the event of intensive use, it is possible to replace sockets 13 and shock absorber 11. During use or cleaning, do not use talcum powder. Do not put the product into contact with corrosive or abrasive elements such as fresh or salt water, acids, dust, sand or other substances.

At least once a year, have the product checked by an orthopaedic technician for proper operation and wear. If you detect any unusual changes in the product's functioning or hear it making noises, stop using it immediately and notify an orthopaedic technician.

RECOMMENDATIONS - WARNINGS

These products should not be used for purposes other than those described in these instructions. They should only be used by one person and under the conditions prescribed. Before use, the product must be adapted by an orthopaedic technician.

When disposing of the product and packaging, strictly comply with the refuse disposal regulations of your community.

In knee joints, there is risk of entrapment, use caution.

As the product features metal components, use caution in the case of magnetic resonance.

WARRANTY

These products comply with European Medical Devices Directive 93/42/EEC (Spanish Royal Decree 1591/2009). They have undergone tests in accordance with international standard EN ISO 22523 for prostheses and orthoses and EN ISO 10328 for lower-limb prostheses.

ORLIMAN, S.L.U. guarantees all of its products, provided that they have not been tampered with or altered from their original form. The company does not guarantee products whose characteristics have been altered through any kind of misuse, breakage or defect. If you notice any defect or abnormality, notify the supplier immediately in order for the product to be replaced.

ORLIMAN S. L. U.
C/ Ausias March, 3 - Pol. Ind. La Pobla-L'Eliaana - Apdo. de correos 49
C.P.: 46185 - La Pobla de Valldona Valencia - España (Spain)
Tel.: +34 96 272 57 04 - Fax: +34 96 275 87 00 - Exports Tel.: +34 96 274 23 33
E-mail: orto@orliman.com - Exports email: export@orliman.com - www.orliman.com

WARRANTY

2 years from purchase

Reference	Batch no.	Date of purchase
<input type="text"/>	<input type="text"/>	<input type="text"/>
Signature and stamp of the orthopaedic establishment		For the warranty to be valid, it is necessary to complete these details.
<input type="text"/>		<input type="text"/>

LB-I-6H09

EDITION: 05/2015

ISO 13485

ISO 9001