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Polycentric pneumatic knee joint

INSTRUCTIONS FOR USE AND MAINTENANCE

Dear Customer,
Please read the instructions carefully. Make sure that the user understands the instructions for use, paying particular attention to the safety and maintenance warnings.

APPLICATION

The polycentric pneumatic knee joint is intended for patients weighing up to 100 kg with unilateral or bilateral amputation at any level, including amputation at the knee. The knee joint ensures secure fastening in the stance phase and adaptation to walking speed. It is recommended for patients with second and third group activity levels.

TECHNICAL CHARACTERISTICS

Code	6H23
Material of the main parts	Aluminium alloy
Item weight, kg	Maximum 0.8
For patients weighing	up to 100 kg
Flexion angle	Max. 150°
Joint in the proximal part	Pyramid
Joint in the distal part	Bushing 30 mm
Total height	235 mm
Device height	179 mm
Max. speed	4 Km/h
Warranty period	3 years


DESCRIPTION OF THE STRUCTURE

The knee joint consists of support 1 connected by means of fork 2 and two connections 3 to housing 4. Axes 5 rotate around bearings 6 creating a closed kinematic loop. In a position of maximum stretching, the instantaneous axis of rotation of the assembly is located significantly higher than the knee joint and above the load line, which provides a secure fit for the knee joint in the stance phase.

Housing 4 features a dual-chamber pneumatic cylinder which ensures control of the knee joint in the swing phase. Inside the housing, a spring folding device has been included. Located in the central part of the housing are valves 7 and 8. Valve 7, marked "F", controls flexion resistance and valve 8, marked "E", controls extension resistance.


For connection to other parts of the prosthesis, support 1 features pyramid 9, and the bottom of housing 4 is designed in the shape of a bushing extendible by means of connecting screw 10.

The knee joint is compatible with similar lower-limb prosthetic systems made by other manufacturers. The connection elements comply with the requirements of Russian standard GOST R 51191.

- 1 Support
- 2 Fork
- 3 Connection
- 4 Housing
- 5 Axis
- 6 Bearing
- 7 Flexion adjustment valve
- 8 Extension adjustment valve
- 9 Pyramid
- 10 Connecting screw

FITTING THE KNEE JOINT TO THE PROSTHESIS

RECOMMENDATIONS The knee joint maintains all of its functions when fitted in conjunction with articulated foot 1H 04Y and, for adjusting the length of the leg, use of joint 3H 01A/400, 3H 01C/400 or 3H 01T/400 is recommended.

ATTENTION To correctly assemble the prosthesis, follow the fitting diagram (see Figure 2). Not strictly following the prosthesis fitting diagram could lead to loss of stability in the stance phase, causing the patient to fall over. The assembly of the prosthesis should be performed in two stages. In the first, the prosthesis is adjusted and, in the second, it is fitted directly onto the patient. For the first stage, a laser level is required.

See Figure 2a for how the knee joint should be fitted to the prosthesis. The adjustment axis of the prosthesis on the sagittal plane must pass through the centre of the proximal area of the side of the receiver socket, through the anterior axis of rotation 1 of the knee joint, at a distance of 30 mm from the midline of the foot towards the heel. This fitting diagram should be used for all types of foot joint, regardless of the recommended distance shown in the technical documentation of the foot joint.

On the prosthesis, during adjustment, 5 mm should be added to the height of the heel. The outer rotation of the knee joint should be approximately 5 degrees (see Figure 2b).

The angle between the line that connects the centres of the proximal and distal parts of the receiver socket with respect to the adjustment axis on the sagittal plane is 3 to 5 degrees, however, the individual characteristics of the patient should be taken into account and the position of the receiver socket corrected. Fig. 2b


Figure 2b

ATTENTION If the position of the receiver socket is not taken into account, the knee joint will excessively incline forwards, causing malfunction and unnecessary wear and tear to the knee joint.

RECOMMENDATIONS To correct the position of the knee joint and the receiver socket in the joint, use of adjustable adapter 4H 16T, 4H 16C or 4H 24T, 4H 24C is recommended.


If the amputation is at knee level, it is necessary to use adapter 4H 13T or 4H 13S.

In the second stage, the final adjustment of the prosthesis is carried out in a static position on the same device to balance the weight with the laser level during the test walk.

The patient is positioned in the device to balance the weight. The height of both legs is made equal and the load between them is distributed as evenly as possible. At the same time, the beam of the laser level is extended along the load line (see Figure 2c). Adjustment of the prosthesis is governed exclusively by the change in plantar flexion, so that the load axis reflected by the laser beam passes at a distance of 35 mm from the axis of rotation of the knee joint.

The final adjustment of the prosthesis is carried out through a test walk by the patient and adjustment of the pneumatic cylinder. The stride rate and length of the orthopaedic leg and the healthy leg should be the same or similar.

To complete the assembly of the prosthesis, screw 10 should be tightened with a torque of 10 Nm.

ADJUSTING THE PNEUMATIC CYLINDER

Pneumatic control of the swing phase prevents the leg with the prosthesis from separating excessively during flexion and ensures smoother extension without sudden jolts.

NOTE The knee joint is supplied from the factory with a pneumatic cylinder setting optimised and adapted for most patients.

For any modification, follow these instructions:

During the test walk, adjust the leg's flexion damping.

- Rotate valve "F" clockwise to increase resistance and, consequently, reduce the angle of separation of the leg.

- Rotate valve "F" anticlockwise to decrease resistance and, consequently, increase the angle of separation of the leg.

Then, adjust the leg's extension damping.

- Rotate valve "E" clockwise to increase resistance when the prosthesis extends abruptly and in the event of strong jolts during the final phase of extension of the prosthesis.

- Rotate anticlockwise to decrease resistance when walking, when the prosthesis partially extends. Carefully adjust the ratio between the flexion and extension values.

ATTENTION Extension can be adjusted as long as complete extension of the knee joint is not impeded. If it is impeded, the patient could fall over.

MAINTENANCE

During the warranty period, the knee joint requires no maintenance.

ATTENTION Inform your patient of the following:

Do not expose the prosthesis to agents that can cause corrosion to its metal components, such as fresh or salt water, acids or other liquids, or the effects of external environmental factors that could lead to the deterioration of its components, such as dust and sand.

Usage and environmental conditions can affect the functional characteristics of the knee joint. In order to avoid risks to the health of the patient, the knee joint should not be used if significant changes in its functioning are noticed. These significant changes can include partial extension, difficulty walking, problems controlling the swing phase and stability in the stance phase, squeaking noises or jolts, etc.

Solution: contact a specialist technical service to check the functioning of the prosthesis.

WARRANTY

These products comply with European Medical Devices Directive 93/42/EEC (Spanish Royal Decree 1591/2009). They have undergone tests in accordance with international standard EN ISO 22523 for prostheses and orthoses and EN ISO 10328 for lower-limb prostheses.

ORLIMAN, S.L.U. guarantees all of its products, provided that they have not been tampered with or altered from their original form. The company does not guarantee products whose characteristics have been altered through any kind of misuse, breakage or defect. If you notice any defect or abnormality, notify the supplier immediately in order for the product to be replaced.

WARRANTY

4 Km/h from purchase

Reference	Batch no.	Date of purchase
Signature and stamp of the orthopaedic establishment		For the warranty to be valid, it is necessary to complete these details.