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6H24

› Single-axis pneumatic knee joint

INSTRUCTIONS FOR USE AND MAINTENANCE

Dear Customer,
Please read the instructions carefully. Make sure that the user understands the instructions for use, paying particular attention to the safety and maintenance warnings.

APPLICATION

The single-axis pneumatic knee joint is indicated for patients with unilateral or bilateral amputation at any level, except for amputation at the knee. The pneumatic knee joint ensures a secure stance phase and adaptation to walking speed. It is recommended for patients with second and third group activity levels.

TECHNICAL CHARACTERISTICS


Code	6H24
Material of the main parts	Aluminium alloy
Item weight, kg	Maximum 0.95
For patients weighing	up to 125 kg
Load capacity in accordance with current Russian legislation (GOST R 51191)	R6
Flexion angle	Max. 120°
Max. speed	6 Km/h
Warranty period	3 years

DESCRIPTION OF THE STRUCTURE


The single-axis pneumatic knee joint consists of housing 1 (Fig. 1), attached to support 2 by means of fixing piece 3. In movement, towards housing 1 and support 2, pneumatic shock absorber 4 is reinforced, which ensures control of the knee joint in the forward movement phase. The valve's screws are located on the upper cover of the pneumatic shock absorber. Screw 5 regulates flexion shock absorption, screw 6 regulates extension shock absorption and screw 7 regulates additional shock absorption. This is achieved by limiting the angular velocity of the tibial prosthesis from 10° to full extension, which prevents abruptness in extension.

During the stance phase, it is reinforced by the weight of the patient, which is transmitted through base 2 to fixing piece 3. The mechanism of fixing piece 3 prevents rotation. During the forward movement phase, it is unlocked.

The knee joint is compatible with similar lower-limb prosthetic systems made by other manufacturers. The connection elements comply with the requirements of Russian standard GOST R 51191.

- 1 Housing
- 2 Support
- 3 Fixing unit
- 4 Pneumatic shock absorber
- 5 Screw to regulate flexion
- 6 Screw to regulate extension
- 7 Screw to regulate additional extension
- 8 Screw
- 9 Screw
- 10 Axis
- 11 Stop

ADJUSTMENT

The knee joint is supplied with factory settings and is ready for use. It is recommended that the patient starts using it slowly and carefully with a walk to check that it is operating properly. If necessary, the settings can be adjusted as required by the patient.

Adjustment of the knee joint in the stance phase

The sensitivity of fixing unit 3 is adjusted by turning screw 9: turning it clockwise decreases sensitivity; anti-clockwise increases it. To check the safety lock, give screw 9 a ¼ turn and get the patient to try it and indicate a comfortable fixing and unlocking point. The maximum rotation of screw 9 cannot penetrate more than 1.5 mm (or a 1.5 turn) with respect to the plane of fixing piece 3. It should not protrude more than 0.5 mm (or a 0.5 turn).

In the event of intensive use, natural wear occurs in the parts of fixing unit 3 (increased gap) which can cause a decrease in sensitivity due to deterioration of fixing unit 3 in the stance phase. To restore the factory setting, it is necessary to:

- Turn screw 8 clockwise (moving stop 11) until flexion (extension) is rigid.
- Half turn screw 8 anti-clockwise until the rigidity disappears (rubbing).
- Get the patient to take a walk to check that the fixing is correct and the unlocking is comfortable.

Corrective adjustment of the pneumatic shock absorber (forward movement phase)

Problem	Solution	Actions to solve the problem
The extent of flexion is too large	Increase shock absorption	Turn screw 5 clockwise to adjust the flexion
The extent of flexion is limited	Reduce shock absorption	Turn screw 5 anti-clockwise to adjust the flexion
The extension movement is too fast	Increase shock absorption	Turn screw 6 clockwise to adjust the extension
At the moment of maximum extension, the support elements make contact.	Increase shock absorption	Turn screw 7 clockwise to adjust the extension. The optimal difference between the positions of screws 7 and 6 is +2.5 to +3.5 turns.

The extension movement is too limited, resulting in incomplete extension of the knee joint before the stance phase.	Reduce shock absorption	Turn screws 6 and 7 clockwise to prevent contact between the support elements. The optimal difference between the positions of screws 7 and 6 is +2.5 to +3.5 turns.
---	-------------------------	--

NOTE:

Screws 5, 6 and 7 have flat grooves and should be adjusted with the flat screwdriver supplied with the knee joint kit.

For adjustment, screws 5, 6 and 7 should be given an ¼ turn.

Pneumatic shock absorber factory setting:

- Number of turns of the screw anti-clockwise – it is screwed in to the maximum.
- screw 5: 4.5 turns
- screw 6: 5.5 turns
- screw 7: 2.5 turns

CAUTION: Adjustment of the pneumatic shock absorber must ensure complete extension of the knee joint before the stance phase. The screws should not be completely turned clockwise – a minimum of one turn should be held in reserve.


FITTING THE KNEE JOINT TO THE PROSTHESIS

The system for fitting the knee joint is shown in Fig. 2. The recommended values, ΔL and H, are shown in the prosthetic foot's accompanying documentation.

The axis must pass through the midpoint of the top of the reception capsule and through the point located 10 mm from the centre of the rotation axis of the knee joint and at ΔL of the midfoot.

To ensure safe adjustment in the stance phase, as well as agility in the transition to the forward movement phase, it is necessary to adjust the position of the prosthesis on the patient taking into consideration the load axis as shown in Figure 2. Adjustment of the prosthesis with respect to the load axis is carried out exclusively at the expense of the plantar flexion of the prosthesis so that the load axis is at a distance of 35-40 mm from the rotation axis of the knee joint.

Carry out the necessary dynamic optimisation adjustments during the test walk.

MAINTENANCE

During use or cleaning, do not use talcum powder. Do not put the product into contact with corrosive or abrasive elements such as fresh or salt water, acids, dust, sand or other substances.

At least once a year, have the product checked by an orthopaedic technician for proper operation and wear. If you detect any unusual changes in the product's functioning or hear it making noises, stop using it immediately and notify an orthopaedic technician.

RECOMMENDATIONS - WARNINGS

These products should not be used for purposes other than those described in these instructions. They should only be used by one person and under the conditions prescribed.

Before use, the product must be adapted by an orthopaedic technician. When disposing of the product and packaging, strictly comply with the refuse disposal regulations of your community.

In knee joints, there is risk of entrapment, use caution.

As the product features metal components, use caution in the case of magnetic resonance.

WARRANTY

These products comply with European Medical Devices Directive 93/42/EEC (Spanish Royal Decree 1591/2009). They have undergone tests in accordance with international standard EN ISO 22523 for prostheses and orthoses and EN ISO 10328 for lower-limb prostheses.

ORLIMAN, S.L.U. guarantees all of its products, provided that they have not been tampered with or altered from their original form. The company does not guarantee products whose characteristics have been altered through any kind of misuse, breakage or defect. If you notice any defect or abnormality, notify the supplier immediately in order for the product to be replaced.

ORLIMAN S. L. U.
C/ Ausias March, 3 - Pol. Ind. La Pobla-L'Eliaana - Apdo. de correos 49
C.P.: 46185 - La Pobla de Valldona Valencia - España (Spain)
Tel.: +34 96 272 57 04 - Fax: +34 96 275 87 00 - Exports Tel.: +34 96 274 23 33
E-mail: orto@orliman.com - Exports email: export@orliman.com - www.orliman.com

WARRANTY

6 Km/h from purchase

Reference	Batch no.	Date of purchase
Signature and stamp of the orthopaedic establishment	For the warranty to be valid, it is necessary to complete these details.	

LB-I-6H24

EDITION: 09/2015

ISO 13485

ISO 9001