# 7H08A

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# 7H08A

## › Hole lock

### INSTRUCTIONS FOR USE AND MAINTENANCE

Dear Customer,  
Please read the instructions carefully. Make sure that the user understands the instructions for use, paying particular attention to the safety and maintenance warnings.

### APPLICATION

The hole lock is used in tibial and femoral prostheses and is designed to fix the silicone liner to the socket, ensuring secure fastening of the prosthesis to the patient's stump without any additional fastening.


### TECHNICAL CHARACTERISTICS

Code	7H08A
Housing material	Aluminium
Bushing and grooved pin material	Stainless steel
Item weight	185 g
For patients weighing	up to 125 kg
Temperature range for use	From -40 °C to +40 °C
Distance between the two pairs of rows of M6 holes for fastening the DCR	36 mm
Warranty period	2 years

### DESCRIPTION OF THE STRUCTURE

The hole lock consists of housing 1 (Fig. 1), pin 2, bushing 3, extension 8 and the locking mechanism. The locking mechanism consists of mechanism housing 4, gear axis 5, free coupling 9, spring 6 and bushing 7.

To prevent the hole lock from unscrewing from housing 1, it features installation screw 11. The hole lock kit also includes screws 10 for fastening to the socket. Socket 8 has an M10 thread for connection to the other components of the hole lock. The hole lock features four M6 holes for fastening the socket support to other compatible lower-limb prosthetic systems made by other manufacturers.

1 Housing	1 unit	7 Bushing	1 unit
2 Grooved pin	1 unit	8 Extension	1 unit
3 Bushing	1 unit	9 Free coupling	1 unit
4 Locking mechanism housing	1 unit	10 Screw	4 unit
5 Gear axis	1 unit	11 Installation screw	1 unit
6 Spring	1 unit		

### FITTING THE ITEM TO THE PROSTHESIS

The fitting of the hole lock is carried out with the 4 screws 10. Screw tightening torque - 6 Nm.

When fitting the hole lock to the prosthesis, it is possible to shorten extension 8 to the area free of grooves.

### OPERATION

To fit the prosthesis, it is necessary to insert the stump with the silicone liner (and pin 2 fastened to it) into the socket. In this case, fixing is automatic when pin 2 attaches to gear axis 5.

If pin 2 has to be adjusted, it is necessary to turn extension 8 clockwise.

To remove the prosthesis, it is necessary to press extension 8 to the stop, which releases pin 2 from its connection to gear axis 5, and then remove the stump from the liner.

If it is necessary to loosen the locking mechanism, unscrew installation screw 11 with the 2.5 mm Allen key.

### ADJUSTMENT AND MAINTENANCE

During the period of use, no adjustment is necessary. All settings are carried out at the factory.

In order for the lock to function reliably for a long time, it is desirable to prevent it from becoming dirty, particularly clogging by sand or other particles, which causes increased wear.

### RECOMMENDATIONS - WARNINGS

These products should not be used for purposes other than those described in these instructions. They should only be used by one person and under the conditions prescribed.

Before use, the product must be adapted by an orthopaedic technician.

When disposing of the product and packaging, strictly comply with the refuse disposal regulations of your community.

In knee joints, there is risk of entrapment, use caution.

As the product features metal components, use caution in the case of magnetic resonance.

### WARRANTY

These products comply with European Medical Devices Directive 93/42/EEC (Spanish Royal Decree 1591/2009). They have undergone tests in accordance with international standard EN ISO 22523 for prostheses and orthoses and EN ISO 10328 for lower-limb prostheses. ORLIMAN, S.L.U. guarantees all of its products, provided that they have not been tampered with or altered from their original form. The company does not guarantee products whose characteristics have been altered through any kind of misuse, breakage or defect. If you notice any defect or abnormality, notify the supplier immediately in order for the product to be replaced.

ORLIMAN S. L. U.  
C/ Ausias March, 3 - Pol. Ind. La Pobla-L'Eliaana - Apdo. de correos 49  
C.P.: 46185 - La Pobla de Valbona Valencia - España (Spain)  
Tel.: +34 96 272 57 04 - Fax: +34 96 275 87 00 - Exports Tel.: +34 96 274 23 33  
E-mail: orto@orliman.com · Exports email: export@orliman.com · www.orliman.com

### WARRANTY

2 years from purchase

Reference	Batch no.	Date of purchase
<input type="text"/>	<input type="text"/>	<input type="text"/>
Signature and stamp of the orthopaedic establishment		For the warranty to be valid, it is necessary to complete these details.
<input type="text"/>		<input type="text"/>

LB-1-7H08A

EDITION: 05/2015

ISO 13485

ISO 9001