8H04CA - 8H04CA/65 - 8H04C - 8H04CT

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ISO 9001

8H04CA - 8H04CA/65 - 8H04C - 8H04CT

3-prong rotatable female socket adapter

INSTRUCTIONS FOR USE AND MAINTENANCE

Dear Customer,
Please read the instructions carefully. Make sure that the user understands the instructions for use, paying particular attention to the safety and maintenance warnings.

APPLICATION

The 3-prong rotatable female socket adapter is used in femoral or tibial prostheses. The 3-prong rotatable female socket adapter is fitted to the socket and serves to connect it to other modules of the prosthesis.

TECHNICAL CHARACTERISTICS


Code	8H04CA	8H04CA/65	8H04C	8H04CT
Material	Stainless steel			
	Aluminium	Stainless steel	Titanium	
Item weight	157 g	285 g	204 g	173 g
Radius of the spherical section, R	50 mm			
Temperature range for use	From -40 °C to +40 °C			
Enables angular adjustment (of the whole prosthesis)	on the frontal and sagittal planes in a range of not less than $\pm 7.5^\circ$			
	on the horizontal plane of not less than $\pm 15^\circ$			
Warranty period	3 years			

DESCRIPTION OF THE STRUCTURE

The 3-prong-rotatable female socket adapter consists of anchor 1 (see Fig. 1) and adjustment pyramid 2, which is screwed to the first by means of the M3x1.5 thread. Anchor 1 features a slot which is tightened with screw 3 to fasten adjustment pyramid 2. On the prosthesis, this connection enables adjustment of the rotation angle on the horizontal plane in a range of not less than $\pm 15^\circ$.

Adjustment pyramid 2 features female spherical surface bearing A and is supplied with two pairs of screws 4 and 5 of 12 mm and 16 mm in length for connection of the socket to the other modules of the prosthesis. This connection enables angular adjustment on the frontal and sagittal planes in a range of not less than $\pm 7.5^\circ$.

Anchor 1 has been manufactured by casting followed by subsequent machining of the piece. The 3-prong rotatable female socket adapter is compatible with Orliman's modular lower limb prosthetic systems.

1 Anchor	1 unit	4 M8 screw, L=12 mm	2 units
2 Adjustment pyramid	1 unit	5 M8 screw, L=16 mm	2 units
3 Screw	1 unit		

FITTING THE ITEM TO THE PROSTHESIS

The 3-prong rotatable female socket adapter is fitted to the liner by means of lamination.


ATTENTION: Screw 3 (Fig. 1) must be located on the sagittal plane, towards the outer side of the prosthesis. Incorrect fitting of the anchor causes damage to the joint between the adjustment pyramid and the anchor.

During assembly of the 3-prong rotatable female socket adapter, it is possible to flex the prongs of the anchor to a radius of not less than R15 mm.

When fitting adjustment pyramid 2 to anchor 1 (and when performing the rotation during correction of the construction schema of the prosthesis), it is necessary to ensure that space S (Fig. 1) between surface B of adjustment pyramid 2 and the end of anchor 1 is not more than 1 mm ($S < 1$ mm).

When connecting adjustment pyramid 2 to other modules of the prosthesis, it is necessary to place the short screws where the space between the ring and the pyramid is smaller, and the long screws, where the space is greater, see Fig. 2.

During final assembly of the prosthesis, it is recommended to apply Loctite 241 threadlocker to screws 4 and 5 and tighten them with a torque of 15 Nm and screw 3 with a torque of 8 Nm.

RECOMMENDATIONS - WARNINGS

These products should not be used for purposes other than those described in these instructions. They should only be used by one person and under the conditions prescribed.

Before use, the product must be adapted by an orthopaedic technician. When disposing of the product and packaging, strictly comply with the refuse disposal regulations of your community.

In knee joints, there is risk of entrapment, use caution. As the product features metal components, use caution in the case of magnetic resonance.

WARRANTY

These products comply with European Medical Devices Directive 93/42/EEC (Spanish Royal Decree 1591/2009). They have undergone tests in accordance with international standard EN ISO 22523 for prostheses and orthoses and EN ISO 10328 for lower-limb prostheses. ORLIMAN, S.L.U. guarantees all of its products, provided that they have not been tampered with or altered from their original form. The company does not guarantee products whose characteristics have been altered through any kind of misuse, breakage or defect. If you notice any defect or abnormality, notify the supplier immediately in order for the product to be replaced.

ORLIMAN S. L. U.
C/ Ausias March, 3 - Pol. Ind. La Pobla-L'Eliaana - Apdo. de correos 49
C.P.: 46185 - La Pobla de Vallbona Valencia - España (Spain)
Tel.: +34 96 272 57 04 - Fax: +34 96 275 87 00 - Exports Tel.: +34 96 274 23 33
E-mail: orto@orliman.com · Exports email: export@orliman.com · www.orliman.com

WARRANTY

3 years from purchase

Reference	Batch no.	Date of purchase
<input type="text"/>	<input type="text"/>	<input type="text"/>
Signature and stamp of the orthopaedic establishment		For the warranty to be valid, it is necessary to complete these details.
<input type="text"/>		<input type="text"/>

LB-I-8H04

EDITION: 05/2015

ISO 13485

ISO 9001