8H05AC - 8H05C - 8H05CT

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ISO 9001

8H05AC - 8H05C - 8H05CT

3-prong rotatable male socket adapter

INSTRUCTIONS FOR USE AND MAINTENANCE

Dear Customer,
Please read the instructions carefully. Make sure that the user understands the instructions for use, paying particular attention to the safety and maintenance warnings.

APPLICATION

The 3-prong rotatable male socket adapter is used in tibial or femoral prostheses. The 3-prong rotatable male socket adapter is fitted to the socket and serves to connect it to other modules of the prosthesis.

TECHNICAL CHARACTERISTICS

Code		8H05AC	8H05C	8H05CT
Material	anchor	Stainless steel		
	adjustment pyramid	Aluminium+ Stainless steel	Stainless steel	Titanium
Item weight		174 g	199 g	162 g
For patients weighing		up to 100 kg		
Radius of the spherical section, R		50 mm		
Temperature range for use		From -40°C to +40°C		
Enables angular adjustment (of the whole prosthesis)		on the frontal and sagittal planes in a range of not less than $\pm 7.5^\circ$		
		on the horizontal plane of not less than $\pm 15^\circ$		
Warranty period		3 years		

DESCRIPTION OF THE STRUCTURE

The 3-prong rotatable male socket adapter consists of anchor 1 (see Fig. 1) and adjustment pyramid 2, which is screwed to the first by means of the M36x1.5 thread. Anchor 1 features a slot which is tightened with screw 3 to fasten adjustment pyramid 2. On the prosthesis, this connection enables adjustment of the rotation angle on the horizontal plane in a range of not less than $\pm 15^\circ$.

Adjustment pyramid 2 features spherical surface bearing A and pyramid B for connection to other modules of the prosthesis. This connection enables angular adjustment on the frontal and sagittal planes in a range of not less than $\pm 7.5^\circ$.

Anchor 1 is manufactured by casting followed by machining of the piece.

The 3-prong rotatable male socket adapter includes a protective cap designed to prevent the resin from making contact with the threaded surface of anchor 1 during lamination of the socket and also includes an adjustable nylon core with a screw to fit it to the anchor.

The 3-prong rotatable male socket adapter is compatible with Orliman's modular lower-limb prosthetic systems.

FITTING THE ITEM TO THE PROSTHESIS

The 3-prong rotatable male socket adapter is fitted to the socket by lamination. In addition, the socket mould is firstly fastened to the expandable adjustment pyramid with a screw and then the anchor, in which the threaded hole is closed with a cover, is fitted.

ATTENTION: Screw 3 (Fig. 1) must be located on the sagittal plane, towards the outer side of the prosthesis.

Incorrect fitting of the anchor causes damage to the joint between the adjustment pyramid and the anchor.

During assembly of the 3-prong rotatable male socket adapter, it is possible to flex the prongs of the anchor to a radius of not less than R15 mm.

When fitting adjustment pyramid 2 to anchor 1 (and when performing the rotation during correction of the construction schema of the prosthesis), it is necessary to ensure that space S (Fig. 1) between surface C of adjustment pyramid 2 and the end of anchor 1 is not more than 1 mm (S<1 mm).

During final assembly of the prosthesis, screw 3 must be tightened with a torque of 8 Nm.

RECOMMENDATIONS - WARNINGS

These products should not be used for purposes other than those described in these instructions. They should only be used by one person and under the conditions prescribed.

Before use, the product must be adapted by an orthopaedic technician.

When disposing of the product and packaging, strictly comply with the refuse disposal regulations of your community.


In knee joints, there is risk of entrapment, use caution.

As the product features metal components, use caution in the case of magnetic resonance.

WARRANTY

These products comply with European Medical Devices Directive 93/42/EEC (Spanish Royal Decree 1591/2009). They have undergone tests in accordance with international standard EN ISO 22523 for prostheses and orthoses and EN ISO 10328 for lower-limb prostheses.

ORLIMAN, S.L.U. guarantees all of its products, provided that they have not been tampered with or altered from their original form. The company does not guarantee products whose characteristics have been altered through any kind of misuse, breakage or defect. If you notice any defect or abnormality, notify the supplier immediately in order for the product to be replaced.

- | | | | |
|----------------------|--------|---------|--------|
| 1 Anchor | 1 unit | 3 Screw | 1 unit |
| 2 Adjustment pyramid | 1 unit | | |

ORLIMAN S. L. U.
C/ Ausias March, 3 - Pol. Ind. La Pobla-L'Eliaua - Apdo. de correos 49
C.P.: 46185 - La Pobla de Valbona Valencia - España (Spain)
Tel.: +34 96 272 57 04 - Fax: +34 96 275 87 00 - Exports Tel.: +34 96 274 23 33
E-mail: orto@orliman.com · Exports email: export@orliman.com · www.orliman.com

WARRANTY

3 years from purchase

Reference	Batch no.	Date of purchase
<input type="text"/>	<input type="text"/>	<input type="text"/>

Signature and stamp of the orthopaedic establishment

For the warranty to be valid, it is necessary to complete these details.

LB-I-8H05

EDITION: 05/2015

ISO 13485

ISO 9001