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English

INSTRUCTIONS FOR USE AND PRESERVATION

Dear customer,

Thank you for your trust in this Orliman product. You have acquired a quality product with a high medicinal value. Please, read the instructions carefully. If you have any doubt, get in contact with your doctor, your orthopaedist specialist or our customer service department. Orliman appreciates your choice and wishes you a speedy recovery.

REGULATION

These products satisfy the European Medical Device Directive 93/42/CEE (RD. 1591/2009). A Risk Analysis has been carried out (UNE EN ISO 14971) minimising all the existing risks. The trials have been carried out in accordance with the European regulation UNE-EN ISO 22523 of Prosthesis and Orthosis.

SUITABLE FOR

- Treatment of clavicular fractures, as a figure-of-eight bandage.
- Description of posture. As a method of rest, limiting mobility and preventing gravity-induced capsulo-ligamentous strain by obtaining a functionally correct position.
- It acts as a method of support, shifting the weight of the upper arm to the back, thereby reducing the biomechanical load on the respective joints. Additionally, it acts as a means of positioning the shoulder in a certain position to minimize the action of certain muscle groups which need to be relaxed, and aligning the fractured segments of the injury in the collarbone.

FITTING INSTRUCTIONS

THE PRODUCT SHOULD BE ADJUSTED BY YOUR ORTHOPAEDIC TECHNICIAN OR IN ACCORDANCE WITH DOCTOR'S ORDERS.

In order to obtain a greater degree of therapeutic effectiveness in the different pathologies and to prolong the life and utility of the product, it is fundamental that the correct size adequate for each patient or user is chosen (enclosed in the packet is a diagram with all the sizes and the equivalent in centimetres). An excess compression can produce intolerance, therefore we recommend that the compression be regulated to an optimum degree.

To fit observe the following aspects:

IC-30:

The aim of this device is to provide the maximum shoulder retropulsion possible, for which there are 4 adjustment straps.

Situated on the back plate, between D2 and D5, immediately adjust the traction strap in a balanced manner and exert the force necessary to obtain the desired support. In cases where an anatomical revision is required, gradually increase the tension in the straps to the tolerance limit of the patient, until the correct posture is obtained.

To apply treatment to reduce fractures and align the clavicle, tighten the straps until the clavicle is correctly aligned.

When using the device in hemiplegic patients with a scapulohumeral subluxation, the back plate (interscapular) is placed at the C7 and D1 level. This is achieved by tightening as much as possible then slightly releasing the lower straps.

Once the plate is in position, tighten all the straps to the maximum extent possible. This supplies the force required to raise the dislocated shoulder and reduce the subluxation.

ET-210:

With the patient in the standing position, centre the back plate by placing the upper straps over the shoulders and the lower straps below, on both sides.

Wrap the lower straps around the stomach below the ribs and fasten them using the velcro in the area above the abdomen.

Then place the upper straps over the shoulders and below the armpit, crossing them across the back and wrap them around the sides in the direction of the abdomen.

Tighten the lower straps as much as necessary to get the desired support for the shoulders, then fasten the straps with the Velcro in the area below the abdomen.

PRECAUTIONS

The constructive material is inflammable.

Do not expose these products to situations where they could ignite. If this were the case remove them quickly and use the necessary measures to put it out. The material used is hypoallergenic, however we cannot guarantee 100% that in some determined cases it will not produce allergies, if this occurs, remove the product and see the doctor that prescribed it.

In case of minor irritations due to sweat, we recommend the use of a cotton sock, to separate the skin from contact with the material. For small inconvenience caused by sweat, we recommend using an interface to separate the skin from contact with the tissue. If rash, irritation or swelling to remove the product and consult a doctor or prosthetist. Contraindicated in open scars with swelling, redness and accumulation heat. The product contains natural rubber latex and may elicit allergic responses in individuals who are sensitised to latex.

RECOMMENDATION-WARNINGS

These products are to be used only by those which have been medically prescribed to do so. Persons unconnected to such prescriptions should not use them. When the product is not being used, store in its original packaging. In order to dispose of the packaging and the product, strictly comply with the legal norms of your community.

The use of the orthosis is conditioned to the recommendations of the prescribing doctor, therefore it must not be used for any other means than those prescribed.

For the validity of the guarantee, the issuing establishment should comply with this sheet of instructions.

MANUFACTURE-CHARACTERISTICS

The material features used in production have been tested and approved and meet all the European quality requirements.

All the products are produced with high quality materials and offer an unbeatable comfort and quality of use. All products offer restraint, stability and compression for the optimum treatment of the pathologies for which they have been designed.

RECOMMENDATIONS FOR PRESERVATION AND CLEANING

Stick the velcros to each other (if the orthosis has any), wash periodically by hand with a damp cloth and a neutral soap. To dry the orthosis, use a dry towel to absorb most of the dampness and leave to dry at room temperature. Do not hang up or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sun light etc. During use or during washing do not use alcohols, ointments or dissolvent liquids. If the orthosis is not properly dried any detergent residues could irritate the skin and deteriorate the product.

GUARANTEES

ORLIMAN, S.L.U. Guarantees all its products, as long as these have not been manipulated or altered from their original state. Does not guarantee those products in which as a result of misuse, deficiencies or breakage of any kind, its original state has been altered. If you observe any deficiency or anomaly, communicate this immediately to the establishment from which it was obtained in order for it to be changed.

ORLIMAN, S.L.U. Appreciates your choice and wishes you a speedy recovery.

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