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English**INSTRUCTIONS FOR USE AND PRESERVATION**

Dear customer,

Thank you for your trust in this Orliman product. You have acquired a quality product with a high medicinal value. Please, read the instructions carefully. If you have any doubt, get in contact with your doctor, your orthopaedist specialist or our customer service department. Orliman appreciates your choice and wishes you a speedy recovery.

REGULATION

These products satisfy the European Medical Device Directive 93/42/CEE (RD. 1591/2009). A Risk Analysis has been carried out (UNE EN ISO 14971) minimising all the existing risks. The trials have been carried out in accordance with the European regulation UNE-EN ISO 22523 of Prosthesis and Orthosis.

SUITABLE FOR

Recommended for patients with tibial amputations.

FITTING INSTRUCTIONS

THE PRODUCT SHOULD BE ADAPTED BY YOUR ORTHOPAEDIC TECHNICIAN OR DOCTOR.

To obtain the greatest degree of therapeutic effectiveness and prolong the useful life of the product, it is fundamental that every patient or user makes the right choice of the most suitable size (a table of sizes is enclosed for guidance in the packing).

To fit observe the following aspects:

Selecting the size (A)

Selecting the size is of utmost importance because a wrong size could cause perspiration problems and a lack of precision in adjusting the prosthesis and while moving; the stabilising matrix and the distal umbrella has to be fitted to the stump exactly. This is why Orliman has a wide range of sizes, in this way covering the needs of the different morphologies of the tibial stumps (see table of sizes).

To select the right size take the measurement at the stump of the patient 4 cm away from the distal end, and then select the size below the measurement obtained, carrying out a test on the patient so as to do a check prior to the casting process. In the case of there being an amputee who is having his first prosthetic device fitted, the swift reduction in the girth of his stump may make it necessary to select two sizes below the liner size.

(Never choose a size above the measurement or round it upwards).

Example: Circumference measurement of the stump at 4cm: 24, you should choose the 23.5 cm size).

Fitting (B)

To carry out the fitting of the liner on the stump, turn it completely inside out as if it were a sock, keeping the textile surface on the inside and medical silicone on the outside, checking that it is clean and dry. Then place the distal cup of the liner on the distal end of the stump, unrolling it over the stump and making sure there are no wrinkles or air pockets and that the fabrics are properly in place. The distal umbrella has to be as close to the centre as possible, enabling the fastening pin to be properly aligned. Then fit the casing in place using the shuttle lock or any other system of suspension that is selected, in this way obtaining the proper and secure adaptation to the prosthesis. Our umbrella system make it possible to adapt any standard pin on the market for use as it is compatible with all of them.

Avoid tugging on the upper end while it is being fitted, preventing tensions on the skin and soft tissue, and also avoid using sharp objects or digging your nails in which could damage the liner. Remember that it is unnecessary to use talcum powder, spray or any other additional substance in fitting this.

PRECAUTIONS

The construction material is inflammable.

Do not expose the products to situations that could cause them to ignite. Should this happen, quickly remove them and use the proper methods for putting this out. The material used is hypoallergenic; however, it is not possible to guarantee completely that it will not cause allergies in certain cases. In the event that it does, remove the product and consult the prescribing doctor.

RECOMMENDATION-WARNINGS

The liner must be carefully washed 2-3 times before being used for the first time.

The pin must be meticulously and securely attached at the fastening point. The cut of the liner must not display sharp edges, angles or folds, because the user could suffer cutaneous reactions and/or cracks may appear in the material.

Protect against direct sources of heat.

The product must only be used by the individual who features in the medical prescription and it is conditional upon the recommendations of the prescribing doctor. It should not be used by anyone to whom it has not been prescribed. When the product is not in use, keep it in its original packaging. When disposing of the product and packaging, strictly comply with the refuse disposal regulations of your community.

In order for the guarantee to be valid, the supplier must fill in the instructions sheet.

MANUFACTURE-CHARACTERISTICS

The materials use in the manufacturing have been tested and approved in all of their characteristics and they meet the requirements and European Regulations on quality that are laid down. Orliman liners are manufactured in medical silicone and are easy to clean.

RECOMMENDATIONS FOR PRESERVATION AND CLEANING

The useful life of the liner will largely depend on the care that we take and the hygiene that we provide to it, as well as properly maintaining it. This is the direct responsibility of the user, who must strictly follow the manufacturer's instructions.

It is recommended that the liner is washed every day using a neutral soap. To facilitate this process, turn the liner inside out, wash the silicon surface with soap, rinse it thoroughly with cold water and dry it with a cloth that is free of fluff, and in the same way maintain continuous hygiene on the stump every time the prosthesis is taken off.

Remember that skin hygiene is very important, and so the skin has to be kept in the best condition, clean and well-hydrated.

Many common bathing and care products, such as soaps, deodorants, colognes, sprays, etc. may cause irritations to the skin.

Check the condition of your liner periodically and report any incident to your technician.

The liner must NEVER be washed in a washing machine.

GUARANTEES

ORLIMAN, S.L.U. Guarantees all its products, as long as these have not been manipulated or altered from their original state. Does not guarantee those products in which as a result of misuse, deficiencies or breakage of any kind, its original state has been altered. If you observe any deficiency or anomaly, communicate this immediately to the establishment from which it was obtained in order for it to be changed.

ORLIMAN, S.L.U. Appreciates your choice and wishes you a speedy recovery.

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