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## English

**INSTRUCTIONS FOR USE AND PRESERVATION**

Dear customer,

Thank you for your trust in this Orliman product. You have acquired a quality product with a high medicinal value. Please, read the instructions carefully. If you have any doubt, get in contact with your doctor, your orthopaedist specialist or our customer service department. Orliman appreciates your choice and wishes you a speedy recovery.

**REGULATION**

These products satisfy the European Medical Device Directive 93/42/CEE (RD. 1591/2009). A Risk Analysis has been carried out (UNE EN ISO 14971) minimising all the existing risks. The trials have been carried out in accordance with the European regulation UNE-EN ISO 22523 of Prosthesis and Orthosis.

**SUITABLE FOR**

Ref.: SD100 STAR BRACE DYNAMIC SHORT SEMI-RIGID BACK SUPPORT.

Ref.: SD101 STAR BRACE DYNAMIC HIGH SEMI-RIGID BACK SUPPORT.

Recommended uses:

Lower back pain, lumbosciatic pain, degenerative processes, muscular debility and weakness, pre- and post-surgery, disc hernias, spondyloarthrosis, strong sacrolumbar support, perfect anatomical adjustment and comfort in day-to-day use.

Ref.: SD103 STAR BRACE DYNAMIC FIX SHORT RIGID CORSET.

Ref.: SD104 STAR BRACE DYNAMIC FIX HIGH RIGID CORSET.

Recommended uses:

Disc hernias. Osteoporosis. Spondyloarthrosis. Spondylolysis with or without stability. Traumatism. Compression fractures. Degenerative disc disease. Post-operation fusion. Chronic muscular disease. Vertebral metastasis treatment.

**FITTING INSTRUCTIONS**

**THE PRODUCT SHOULD BE ADJUSTED BY YOUR ORTHOPAEDIC TECHNICIAN OR IN ACCORDANCE WITH DOCTOR'S ORDERS.**

In order to obtain a greater degree of therapeutic effectiveness in the different pathologies and to prolong the life and utility of the product, it is fundamental that the correct size adequate for each patient or user is chosen (enclosed in the packet is a diagram with all the sizes and the equivalent in centimetres). An excess compression can produce intolerance, therefore we recommend that the compression be regulated to an optimum degree.

**To fit observe the following aspects:**

**Adjusting and assembling the thermoplastic module and/or lumbar cushion:**

- (A-B) The Star Brace Dynamic corset is fitted with a thermoplastic lumbar cushion. Slide the corset using the guide located on the lumbar plate for removal or positioning.
- (C-D) The Star Brace Dynamic Fix corset is fitted with a thermoplastic lumbar module.

**If the thermoplastic module needs adjusting:**

- Remove the lumbar cushion if fitted. If it is not fitted, remove the module by sliding along the fabric cover guides.
- Check the thermoplastic module on the patient and adjust and modify using the hot air gun. Then assemble the thermoplastic module on the fabric cover by introducing the sacrolumbar corset along the rear guides located between the thermoplastic module and the cover, ensuring it is as aligned as possible, and fasten using the Velcro straps.
- Coinciding with the abdominal zone of the corset lapel is a black elastic band which can be used to secure the thermoplastic abdominal plate, protected by a padded fabric cover. Introduce this strap in the groove fitted for this purpose in the thermoplastic piece, ensuring it is properly in place.
- Once all the components have been assembled, adjust it on the patient.

**Adjusting the StarBrace Dynamics (SD100/101) to the patient:**

Place the corset on the patient, ensuring the lumbar cushion is aligned, then fasten the corset and adjust using the Velcro fasteners on the front.

Finally, adjust the required compression by tightening the left and right cord sections, which will pull in the upper and lower parts of the orthosis, and secure using the colour-coded Velcro tabs.

**Adjusting the StarBrace Dynamics Fix (SD103/104) to the patient:**

Once the thermoplastic module fitted with the corset has been adjusted and/or modified, place the corset on the patient, ensuring that the thermoplastic lumbar support is correctly positioned, then wrap the fabric around the patient and adjust using the Velcro fasteners on the front, making sure the abdominal thermoplastic plate is perfectly aligned.

Finally, adjust the required compression by tightening the left and right cord sections, which will pull in the upper and lower parts of the orthosis, and secure using the colour-coded Velcro tabs.

**PRECAUTIONS**

The constructive material is inflammable.

Do not expose these products to situations where they could ignite. If this were the case remove them quickly and use the necessary measures to put it out. The material used is hypoallergenic, however we cannot guarantee 100% that in some determined cases it will not produce allergies, if this occurs, remove the product and see the doctor that prescribed it.

In case of minor irritations due to sweat, we recommend the use of a cotton sock, to separate the skin from contact with the material. For small inconvenience caused by sweat, we recommend using an interface to separate the skin from contact with the tissue. If rash, irritation or swelling to remove the product and consult a doctor or prosthetist. Contraindicated in open scars with swelling, redness and accumulation heat.

**RECOMMENDATION-WARNINGS**

These products are to be used only by those which have been medically prescribed to do so. Persons unconnected to such prescriptions should not use them. When the product is not being used, store in its original packaging. In order to dispose of the packaging and the product, strictly comply with the legal norms of your community.

The use of the orthosis is conditioned to the recommendations of the prescribing doctor, therefore it must not be used for any other means than those prescribed.

For the validity of the guarantee, the issuing establishment should comply with this sheet of instructions. Includes metal components so extreme caution to other uses, and in case of Magnetic Resonance.

**MANUFACTURE-CHARACTERISTICS**

The material features used in production have been tested and approved and meet all the European quality requirements.

All the products are produced with high quality materials and offer an unbeatable comfort and quality of use. All products offer restraint, stability and compression for the optimum treatment of the pathologies for which they have been designed.

**RECOMMENDATIONS FOR PRESERVATION AND CLEANING**

Stick the velcros to each other (if the orthosis has any), wash periodically by hand with a damp cloth and a neutral soap. To dry the orthosis, use a dry towel to absorb most of the dampness and leave to dry at room temperature. Do not hang up or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sun light etc. During use or during washing do not use alcohols, ointments or dissolvent liquids. If the orthosis is not properly dried any detergent residues could irritate the skin and deteriorate the product.

**GUARANTEES**

ORLIMAN, S.L.U. Guarantees all its products, as long as these have not been manipulated or altered from their original state. Does not guarantee those products in which as a result of misuse, deficiencies or breakage of any kind, its original state has been altered. If you observe any deficiency or anomaly, communicate this immediately to the establishment from which it was obtained in order for it to be changed.

ORLIMAN, S.L.U. Appreciates your choice and wishes you a speedy recovery.

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