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English**INSTRUCTIONS FOR USE AND PRESERVATION**

Dear customer,

Thank you for your trust in this Orliman product. You have acquired a quality product with a high medicinal value. Please, read the instructions carefully. If you have any doubt, get in contact with your doctor, your orthopaedist specialist or our customer service department. Orliman appreciates your choice and wishes you a speedy recovery.

REGULATION

These products satisfy the European Medical Device Directive 93/42/CEE (RD. 1591/2009). A Risk Analysis has been carried out (UNE EN ISO 14971) minimising all the existing risks. The trials have been carried out in accordance with the European regulation UNE-EN ISO 22523 of Prosthesis and Orthosis.

SUITABLE FOR

- Stable osteoporotic fracture.
- Osteoporosis.
- Muscle weakness.
- Scheuermann's juvenile kyphosis.
- Kyphosis with back pain.
- Vertebral misalignments.
- Postural rehabilitation.

FITTING INSTRUCTIONS

THE PRODUCT SHOULD BE ADJUSTED BY YOUR ORTHOPAEDIC TECHNICIAN OR IN ACCORDANCE WITH DOCTOR'S ORDERS.

An excess compression can produce intolerance, therefore we recommend that the compression be regulated to an optimum degree.

To fit observe the following aspects:

Wearing lightweight clothing underneath the brace can improve comfort.

A- Fitting the brace to the patient:**A.1-Adjusting the metal plate.**

- 2-Peel off the micro hook closure of the lower textile component of the rear area.
- 3-Separate both components, removing the metal plate inside the lower component.
- 4-Remove the closure in a micro-hook that fixes the plate to the upper textile component to access the complete metal plate.
- 5-With the patient standing and looking forwards, position the stabiliser plate on the patient's back so that the upper edge is approximately 4 cm below C7 and the lower edge is at the level of the coccyx.
- 6-If it is necessary to adjust the height of the plate, remove the 4 screws with a 2.5 mm Allen key, slide both pieces to the desired height and re-fix all of the screws to the ends of the upper and lower plates.
- 7-Adjust the plate to the anatomical shape of the patient's spine beginning with the lumbar lordosis. It is not necessary to use tools for this process.
- 8-Once the adaptation of the metal plate has been completed, you must insert it into the upper textile component and close the respective micro-hook.
- 9-Insert the lower part of the plate into the lower textile component and close the respective micro-hook.

A.1-Adjusting the straps.

- 2-After closing the brace in the abdominal area, begin by adjusting the tension and length of the 4 straps that connect this part to the rear plate. Special attention should be paid to the lower straps, which should be tight. To adjust the length of the straps, release the Velcro fasteners fixed to their ends, trim the straps to the desired length and reattach the fasteners.
- 3-With the patient standing and looking forwards, adjust the tension and length of the 2 shoulder straps in order obtain a correct posture when they are fixed to the abdominal area.
- 4-The intermediate actuation point of these straps can be adjusted in the rear area by moving the Velcro piece that joins the two buckles.

A.1-Adjusting the armpit protectors.

- 2-To move the armpit protectors, detach the 2 Velcro fasteners on their ends and move the protectors until they are in contact with the front area of the shoulder and the armpit area to prevent discomfort from rubbing.
- 3-After situating them correctly, fix their position with the Velcro fasteners on the ends.

A- Placement of the brace by the patient:

After the orthopaedic technician has carried out the fitting, it will not be necessary to move the positions of the straps. For the placement of the brace, only the abdominal closure needs to be handled.

- 1-With the abdominal closure open, insert the arms through the openings of the shoulder straps, as if it were a backpack.
- 2-Insert the fingers into the respective loops on the left and right-hand sides of the closure.
- 3-Press the left-hand side of the closure on the abdomen and pull the right-hand side until the Velcro can be closed.
- 4-To close the Velcro, firstly remove the fingers from the left-hand loop and then from the right-hand loop.

PRECAUTIONS

The constructive material is inflammable.

Do not expose these products to situations where they could ignite. If this were the case remove them quickly and use the necessary measures to put it out. The material used is hypoallergenic, however we cannot guarantee 100% that in some determined cases it will not produce allergies, if this occurs, remove the product and see the doctor that prescribed it.

In case of minor irritations due to sweat, we recommend the use of a cotton sock, to separate the skin from contact with the material. For small inconvenience caused by sweat, we recommend using an interface to separate the skin from contact with the tissue. If rash, irritation or swelling to remove the product and consult a doctor or prosthetist. Contraindicated in open scars with swelling, redness and accumulation heat.

RECOMMENDATION-WARNINGS

These products are to be used only by those which have been medically prescribed to do so. Persons unconnected to such prescriptions should not use them. When the product is not being used, store in its original packaging. In order to dispose of the packaging and the product, strictly comply with the legal norms of your community.

The use of the orthosis is conditioned to the recommendations of the prescribing doctor, therefore it must not be used for any other means than those prescribed.

For the validity of the guarantee, the issuing establishment should comply with this sheet of instructions.


MANUFACTURE-CHARACTERISTICS

The material features used in production have been tested and approved and meet all the European quality requirements.

All the products are produced with high quality materials and offer an unbeatable comfort and quality of use. All products offer restraint, stability and compression for the optimum treatment of the pathologies for which they have been designed.

RECOMMENDATIONS FOR PRESERVATION AND CLEANING

Stick the velcros to each other (if the orthosis has any), wash periodically by hand with a damp cloth and a neutral soap. To dry the orthosis, use a dry towel to absorb most of the dampness and leave to dry at room temperature. Do not hang up or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sun light etc. During use or during washing do not use alcohols, ointments or dissolvent liquids. If the orthosis is not properly dried any detergent residues could irritate the skin and deteriorate the product.

**GUARANTEES**

ORLIMAN, S.L.U. Guarantees all its products, as long as these have not been manipulated or altered from their original state. Does not guarantee those products in which as a result of misuse, deficiencies or breakage of any kind, its original state has been altered. If you observe any deficiency or anomaly, communicate this immediately to the establishment from which it was obtained in order for it to be changed.

ORLIMAN, S.L.U. Appreciates your choice and wishes you a speedy recovery.

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