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English**INSTRUCTIONS FOR USE AND PRESERVATION**

Dear customer,

Thank you for your trust in this Orliman product. You have acquired a quality product with a high medicinal value. Please, read the instructions carefully. If you have any doubt, get in contact with your doctor, your orthopaedist specialist or our customer service department. Orliman appreciates your choice and wishes you a speedy recovery.

REGULATION

These products satisfy the European Medical Device Directive 93/42/CEE (RD. 1591/2009). A Risk Analysis has been carried out (UNE EN ISO 14971) minimising all the existing risks. The trials have been carried out in accordance with the European regulation UNE-EN ISO 22523 of Prosthesis and Orthosis.

SUITABLE FOR

Fractures of the humeral diaphysis, whether by direct or indirect trauma, pathological fractures, surgical treatment of fractures. As an aid to the treatment of upper arm fractures associated with the elbow joint.

FITTING INSTRUCTIONS

THE PRODUCT SHOULD BE ADJUSTED BY YOUR ORTHOPAEDIC TECHNICIAN OR IN ACCORDANCE WITH DOCTOR'S ORDERS.

In order to obtain a greater degree of therapeutic effectiveness in the different pathologies and to prolong the life and utility of the product, it is fundamental that the correct size adequate for each patient or user is chosen (enclosed in the packet is a diagram with all the sizes and the equivalent in centimetres). An excess compression can produce intolerance, therefore we recommend that the compression be regulated to an optimum degree.

To fit observe the following aspects:

Because it is a brace made of thermoplastic, the use of an interface between the patient and brace is recommended (included in the package).

Because the arm is injured, it is recommended not to make any sudden movements with the upper arm when fitting the patient with the brace.

1-Completely open all the brace straps, separating the thermoplastic sections, until there is sufficient room to comfortably insert the arm.

2-Position it so that the mechanical axis is aligned with the anatomical axis, then using a screwdriver unscrew the guides and adjust the brace for the arm and forearm to the length of the patient's arm.

3-After checking that the size is correct, fasten the straps beginning with those nearest to the elbow joint. Then do the same with the wrist and shoulder straps.

Next, pull the chest strap across to the other side, adjust and insert the traction pin until the desired tension is achieved. With a permanent marker, mark the excess thermoplastic, with particular emphasis on the underarm and elbow areas, to enable the required bending. Finally, perform the same operation at the height of the wrist so that it is not too tight and to prevent grazes.

After taking the brace off, remove the surplus thermoplastic.

Modifying the elbow joint:

Open the protective cover of the elbow hinge to gain access to the screws.

Using the key provided, undo the screws and put them in the hinge holes provided to control the mobility of the flex-extension length as required.

Modifying the shoulder strap (arm / forearm support):

Once the articulated brace is adjusted to the upper arm, place one end of the arm support strap in the side area/back of the humerus brace, then position it across the back towards the opposite shoulder. Pull it over the shoulder, with the other end of the strap touching the forearm at wrist level. Adjust the strap by lifting the arm to the desired level.

PRECAUTIONS

The constructive material is inflammable.

Do not expose these products to situations where they could ignite. If this were the case remove them quickly and use the necessary measures to put it out. The material used is hypoallergenic, however we cannot guarantee 100% that in some determined cases it will not produce allergies, if this occurs, remove the product and see the doctor that prescribed it.

In case of minor irritations due to sweat, we recommend the use of a cotton sock, to separate the skin from contact with the material. For small inconvenience caused by sweat, we recommend using an interface to separate the skin from contact with the tissue. If rash, irritation or swelling to remove the product and consult a doctor or prosthetist. Contraindicated in open scars with swelling, redness and accumulation heat. The product contains natural rubber latex and may elicit allergic responses in individuals who are sensitised to latex.

RECOMMENDATION-WARNINGS

These products are to be used only by those which have been medically prescribed to do so. Persons unconnected to such prescriptions should not use them. When the product is not being used, store in its original packaging. In order to dispose of the packaging and the product, strictly comply with the legal norms of your community.

The use of the orthosis is conditioned to the recommendations of the prescribing doctor, therefore it must not be used for any other means than those prescribed.

For the validity of the guarantee, the issuing establishment should comply with this sheet of instructions.

MANUFACTURE-CHARACTERISTICS

The material features used in production have been tested and approved and meet all the European quality requirements.

All the products are produced with high quality materials and offer an unbeatable comfort and quality of use. All products offer restraint, stability and compression for the optimum treatment of the pathologies for which they have been designed.

RECOMMENDATIONS FOR PRESERVATION AND CLEANING

Stick the velcros to each other (if the orthosis has any), wash periodically by hand with a damp cloth and a neutral soap. To dry the orthosis, use a dry towel to absorb most of the dampness and leave to dry at room temperature. Do not hang up or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sun light etc. During use or during washing do not use alcohols, ointments or dissolvent liquids. If the orthosis is not properly dried any detergent residues could irritate the skin and deteriorate the product.

**GUARANTEES**

ORLIMAN, S.L.U. Guarantees all its products, as long as these have not been manipulated or altered from their original state. Does not guarantee those products in which as a result of misuse, deficiencies or breakage of any kind, its original state has been altered. If you observe any deficiency or anomaly, communicate this immediately to the establishment from which it was obtained in order for it to be changed.

ORLIMAN, S.L.U. Appreciates your choice and wishes you a speedy recovery.

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