ESPAÑOL

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ENGLISH**INSTRUCTIONS FOR USE, STORAGE AND WARRANTY**

Dear Customer,

Thank you very much for placing your trust in an Orliman product. Please read the instructions carefully. Keep these instructions and the packaging for future reference. If you have any questions or concerns, please contact your doctor, orthopaedic specialist or our customer service department.

ORLIMAN S.L.U. guarantees all its products as long as the original configuration has not been manipulated or altered except for the intended use as described in these instructions. If the products are used in combination with other products, replacement parts or systems, make sure they are compatible and made by Orliman®. It does not guarantee any products with altered characteristics due to improper use, defects or breakage of any kind. To use the warranty, go to the place where you bought the product as they must complete the warranty information table on this instructions sheet. If any serious incidents related to the product occur, notify Orliman S.L.U. and the corresponding competent authority in your country.

Orliman would like to thank you for choosing this product and hopes you a speedy recovery.

REGULATIONS

MD This article is defined as a class I medical device. A Risk Analysis (UNE EN ISO 14971) has been carried out, minimising the existing risks. Tests have been in accordance with European Regulation UNE-EN ISO 22523 on Prosthetics and Orthotics.

INDICATIONS

- Patellar and mediolateral instability.
- Meniscus injuries.
- Effusion and swelling caused by arthrosis and arthritis.
- Fibromyalgia.
- Sprains and contusions.
- Post-surgical and post-trauma irritation.

FITTING INSTRUCTIONS

For best therapeutic results considering different pathologies and to extend the useful life of the product, it is essential to choose the correct size for each patient or user. Excessive compression may be intolerable; adjusting the compression to be firm yet comfortable is recommended.

If the product must be adapted, any such adjustments must be done by an orthopaedic specialist or healthcare professional legally certified to do so who must make sure the end user or person responsible for fitting the product properly understands how it works and should be used.

When fitting the product, you must adhere to the following instructions:

A ADAPTING THE JOINTS:

If the patient morphology so requires, the joints must be moulded to prevent discomfort and guarantee proper function.

1-Remove both joints from the upper pocket. Mould the joints and check them on the patient's leg, making sure the mechanical joint axis coincides with the anatomic joint. Profile the shape until it is in contact with the knee, respecting the morphology of the femoral condyles. Make sure the ends of the joints do not pressure the thigh and calf.

B PUTTING ON THE KNEE BRACE:

1-Loosen the straps, hold the knee brace with the thumb on the inside of the joints and put the knee brace on the leg.

2-Make sure the patellar hole coincides with the kneecap so the patellar donut surrounds the outer edges of the kneecap. The joints must be parallel in the middle and at the side of the knee with the mechanical axis coinciding with the anatomical axis.

3-Adjust the textile straps beginning with the bottom strap.


4-If the stabilisation provided by the strap or joints is not necessary, they can be removed.


⚠ PRECAUTIONS

Before each use, check that all product components are present as per the fitting process. Periodically check the conditions of the product. If you observe any defect or anomaly, immediately report it to the issuing establishment.

This product is made of inflammable material. Do not expose the products to situations that could set them on fire. In the event of a fire, quickly get them off your body and use the proper resources to extinguish the fire.

To avoid minor discomfort caused by sweating, we recommend using some type of cotton fabric to separate the skin from contact with the product material. For discomfort such as chafing, irritation and swelling, remove the product and see a doctor or orthopaedic specialist. The product should only be used on healthy skin. It is not recommended for use over open scars with swelling, redness or hotspots.

Products marked with the  symbol contain natural rubber latex and can cause allergic reactions in people sensitive to latex.


Products marked with the  symbol contain ferromagnetic components and, therefore, extreme precaution must be taken if you undergo an MRI scan or are exposed to radiation associated with diagnostic or therapeutic procedures.

RECOMMENDATIONS-WARNINGS

The use of these products is conditioned by the indications. Although the product is not defined as a single-use device, using it on a single patient only is recommended and only for the intended purposes as described in these instructions or by a healthcare professional. When disposing of the product and its packaging, you must strictly adhere to the legal regulations in your community.

RECOMMENDATIONS FOR STORAGE AND WASHING

When not using the product, store it in the original packaging in a dry place at room temperature. Stick the Velcro to each other (if the orthotic device has them), frequently wash by hand with warm water (30° C max.) and mild soap. To dry the product, use a dry towel to absorb as much moisture as possible and let it dry at room temperature. Do not hang it up or iron the product and do not expose it to direct heat sources such as stoves, dryers, direct sun exposure, etc. When using or cleaning the product, do not use abrasive or corrosive substances, alcohol, ointments or liquid solvents. If not dried off properly, the detergent residue may irritate the skin and cause the product to deteriorate.

Remove the rigid splint before washing the brace. To do so, open the pocket on the inside of the wrist brace. Once washed and completely dry, re-insert the rigid splint on the proper side depending on which hand is affected, as per the fitting instructions.

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