



The pads redistribute the compression of the bandage from the protruding epicondyles towards the fleshy parts around them.

**Ref.: 8301** BEIGE

**Ref.: 9301** BLACK

## ELASTIC ELBOW SUPPORT WITH VISCOLASTIC PADS

### Characteristics:

Developed in breathable elastic woven fabric using flat weaving, it incorporates two anatomically-shaped Viscolastic pads which redistribute the compression on the epicondyle and epitrochlea.

### Indications:

- Epicondyle pathologies (e.g. tennis elbow, golf elbow), osteoarthritis and arthritis.
- States of postoperative and post-traumatic irritation, tendinopathies, epicondylitis, prevention of sports or work-related injuries.

### Colours:

- 8301: ■ Beige.
- 9301: ■ Black.

SIZES	PERIMETER CM	LENGTH CM
1	17-19	21,5
2	19-21	22
3	21-23	22,5
4	23-25	24
5	25-27	24
6	27-29	24
7	29-31	24,5

12CM FROM THE ELBOW

