



Ref.: EP-20A



EPICONDYLITIS ARMBAND

Characteristics:

Made of velour on both sides and reinforced with a core of polyethylene and plastazote padding. The closure strap is adjusted on the forearm and fastens back on itself, making it individually adaptable according to the pressure exerted by its elasticated cuff with hook fastener.

Indications:

- Sports and work injury prevention.
- Epicondylitis (tennis elbow).
- Epitrochleitis (golfer's elbow).

Colour:

- Blue.

SIZES	ONE SIZE	3CM FROM THE ELBOW
-------	----------	------------------------