



Ref.: FP-D74 RIGHT

Ref.: FP-I74 LEFT

Characteristics:

Immobilising splint in thermoplastic which includes a very comfortable breathable interface. The splint has a small micro-velcro fastener which can be adjusted to fit to the required length.

Indications:

- Tendonitis.
- Painful and inflammatory sequels.
- Instability of metacarpophalangeal joints.
- Sports injuries (eg, skier's thumb).
- Irritations of the metacarpophalangeal joints.
- Immobilisation of the trapezo-metacarpal joint. (Rhizarthrosis of the thumb).
- Immobilisation of the metacarpophalangeal joints in the thumb (after injury or surgery).

Colour:

- Beige.

SIZES	PERIMETER CM	LENGTH CM	
1	<17	11	WRIST
2	17-19	12	
3	>19	13	LENGTH

The breathable lining can be ordered separately. Ref: Forro-FP74