

Thermo-med®



Ref.: 4103-A

NEOPRENE KNEE SUPPORT WITH THIGH AND CALF OPENINGS FLEXIBLE LATERAL STABILISERS AND ADJUSTMENT STRAPS

Characteristics:

4,5 mm neoprene knee support with flexible lateral stabilisers and upper and lower adjustment straps. Interchangeable patellar pad. Easy to fit and remove thanks to the thigh and calf openings.

Indications:


- Light contusions.
- Arthrosis and arthritis.
- Post-surgical treatments. Rehabilitation.
- Thermo-compression, relief, protection, stability, centring of the kneecap, adaptation to specific anatomies.

Colours:


- 4103-A: ■ Blue.
- 5103-A: ■ Beige.

SIZES	PERIMETER CM
1	32-34
2	34-36
3	36-38
4	38-41
5	41-44
6	44-49

Perimeter in cm around centre of kneecap



HEIGHT



34 cm

