



Ref.: 4116

NEOPRENE STABILISER KNEE SUPPORT WITH SILICONE KNEEPAD, LATERAL FLEXIBLE STABILISERS AND ADJUSTMENT STRAPS

Characteristics:

4,5 mm neoprene knee support with open kneecap and spiral reinforcement stays for lateral stabilisation. With silicone kneepad and straps with elastic portion for better adjustment. Rear thigh and calf opening with Velcro for better adaptation.

Indications:

- Functional treatment in distensions, instabilities and contusions.
- Arthrosis and arthritis.
- The upper and lower openings lighten excessive pressures.

Colours:

- 4116: ■ Blue.
- 5116: ■ Beige.

SIZES	PERIMETER CM
1	32-34
2	34-36
3	36-38
4	38-41
5	41-44
6	44-49

Perimeter in cm around centre of kneecap

HEIGHT

34 cm