



4119

## Wraparound stabilising knee brace



### + Characteristics

- Made of closed cell neoprene for adequate body heat conservation.
- Wrap-around design for a perfect fit with all types of morphologies.
- With lateral spiral stabilisers.

### + Recommended for

- Arthrosis and arthritis.
- Tendinitis and bursitis.
- Contusions and muscle pain.
- Chondromalacia patellae (runner's knee).
- Kneecap dislocation.
- Non-specific knee pain which improves with heat.
- Support and protection for sports.
- Post-surgical and post-trauma treatment.

### + Effects

- Thermal compression, stabilisation and support (patella).



One Size



Mín: 25 cm  
Máx: 48 cm



27 cm