



**Ref.: 7117D RIGHT**

**Ref.: 7117I LEFT**

**KNEE SUPPORT FOR LATERAL OR MEDIAL PATELLA CONTROL**

**Characteristics:**

Fabric knee support made from three-layer breathable elastic material. Composed of three layers – elastic microfibre fabric, polyurethane foam and terry cloth – which facilitates rapid dispersion of sweat, greater elasticity and a better fit. Features silicone bubble grip to prevent slippage. Features flexible mediolateral stabilisers, semi-circular silicone patella cushion and patella pad, also semi-circular.

**Indications:**

- External or internal hyperpressure syndromes and patellar subluxations.
- Prevention of recurring injuries (e.g. patellar luxation).
- Post-surgical treatment of the lateral or medial section of the retinaculum.

**Colour:**

- Black.

SIZES	PERIMETER CM	
1	32-34	<p>Perimeter in cm around centre of kneecap</p>
2	34-36	
3	36-38	<p>HEIGHT 35 cm</p>
4	38-41	
5	41-44	
6	44-49	

Pad position	Select ref.
Right outer	7117D (Right)
Right inner	7117I (Left)
Left outer	7117I (Left)
Left inner	7117D (Right)