



Ref.: AB14 BLACK

Ref.: AB24 BEIGE



CALF SUPPORT FOR THE BOXIA® DROP FOOT ANKLE BRACE

Characteristics:

The Boxia drop foot ankle brace's supramalleolar support stores a large amount of energy, similar to a pivot (support point), to enable it to transfer the necessary forces to lift a paralysed foot. In cases in which a patient shows intolerance due to excessive force, the Boxia calf support is recommended. It helps to stabilise the orthosis and restrict forces and stresses by transferring them from the ankle to the calf.

Indications:

- Flacid paralysis.

Colours:

- AB14: ■ Black.
- AB24: ■ Beige.

SIZES AB01	SIZES AB14
0	0
1 - 2	1
3	2

Equivalences for chose the correct size.