



Ref.: OSL1306

Ref.: OSL1307 ●

SOFT ANTI-BEDSORE WRIST PROTECTOR



Ref.: OSL1312

Ref.: OSL1313 ●

SOFT ANTI-BEDSORE MITTEN



Ref.: OSL1304

Ref.: OSL1305 ●

SOFT ANTI-BEDSORE ELBOW PROTECTOR

COMMON CHARACTERISTICS



anti-microbial

Characteristics:

Anti-bedsore protector with anti-bacterial properties. Soft fabric with Sanitized treatment.

It does not absorb moisture, keeping the skin always dry and improving heat control.

With flat seams to prevent possible chafing.

Ambidextrous protector, including straps for easy, adjustable fastening.

Effects:

- Prevents the appearance of pressure ulcers or sores by distributing pressure over bony prominences and other vulnerable/sensitive areas.

- Reduces friction and shearing forces.

- Keeps the skin clean and dry, ensuring maximum comfort, hygiene and breathability.

- Controls moisture.

- Promotes and facilitates patient rehabilitation.

- Provides better preventative and therapeutic care.

Indications:

Prevention and auxiliary treatment of pressure sores.

People who spend long periods of time in the same position with parts of the body in contact with the surface of the bed or subjected to constant friction by bedding or rough surfaces.

Protection areas by product:

- Wrist protector: temporal styloid process and situations of venoclysis.

- Mitten: palmar and dorsal areas of the hand; prevents the patient from causing self-injury by limiting gripping ability.

- Elbow protector: olecranon, epicondyle and epitroclea.

Colour:

OSL1306 / OSL1312/ OSL1304: □ White.

OSL1307 / OSL1313 / OSL1305: ● Grey.

Composition:

- Fabric: 90% polyester, 10% acrylic with Sanitized treatment.

- Adjustable closure system: 100% polyamide.

- Straps (if featured): 100% polyester.

- Buckle (if featured) 100% polyacetal (POM).


WRIST

SIZES ONE SIZE



MITTEN

SIZES	WIDTH	INTERIOR LENGTH
1	5-6	24
2	7-9	25,5
3	9-11	27



ELBOW

SIZES ONE SIZE

FOREARM

