

# Boxia® Plus



**Ref.: AB100**

## BOXIA® PLUS

### Characteristics:

Comprised of:

**A** - One supramalleolar part made of micro-perforated breathable fabric that wraps around the ankle. It features a pad in the back that hugs the Achilles tendon and provides maximum comfort and distributes pressure. There are also relief points offering a massage effect on the soft parts. The central anchoring system can be adjusted to the ankle perimeter to guarantee the most adequate position for each patient. An elastic traction band keeps the ankle in a neutral position and connects the supramalleolar part to the "S" hooks or the hook to the footwear.

The anti-slip hook attached to the footwear and a strap that is secured to the laces prevents undesirable movement.

4 pairs of "S" hooks for greater stability and varum-valgum foot control.

**B** - 1 An elastic traction band keeps the ankle in a neutral position and connects the supramalleolar part to the "S" hooks or the hook to the footwear.

**C** - The anti-slip hook attached to the footwear and a strap that is secured to the laces prevents undesirable movement.

**D** - 4 pairs of "S" hooks for greater stability and varum-valgum foot control.

### Indications:

- Weakness or inability to dorsiflex the ankle (clubfoot) as a result of neurological injuries, congenital alterations or muscular disorders.
- Flaccid paralysis associated with hemiparesis.

### Benefits:

- Keeps the ankle in a functional position to prevent foot collapse.
- Enables a more natural and efficient walking pattern.

### Colour:

AB01:  Black.

SIZES	PERIMETER CM
1	14-20
2	20-26
3	26-32

