

Cervical collars



Ref.: CC2400



SEMI-RIGID COLLAR WITH CHIN SUPPORT (adjustable)

Characteristics:

Made from two overlapping polyethylene pieces, height-adjustable by means of a Velcro binding, the upper piece is designed to accommodate the chin and distribute pressure over its entire surface for optimum comfort. Padded on its upper and lower edges, lined with imitation leather, rear Velcro closure system.

Indications:

- Minor trauma of bony (spinous and transverse processes) and soft parts (whiplash).
- Postoperative care, torticollis, cervicalgia, cervicobrachialgia, compression syndromes, rheumatic and/or degenerative disorders.

Colour:

- White.

SIZES	PERIMETER CM
1	30-35
2	35-40
3	40-45

CIRCUMFERENCE

HIGH

HEIGHT	Front	Min.15 cm - Max. 20 cm
	Rear	Min. 7 cm - Max. 12 cm