

Epitec® Fix



Ref.: EP-21G



EPITEC® FIX EPICONDYLITIS ARMBAND

Characteristics:

The external side is made with highly resistant cushioned materials, internally with a moulded and shaped thermoplastic soul, and with silicone points that stick to by the skin's pressure. Easily regulated with the closing band and the elastic band.

Indications:

- Prevention of sporting and work-related injuries.
- Epicondylitis (tennis elbow).
- Epitrocleitis (golfer's elbow).


Colour:

- Black.

SIZES

ONE SIZE

3CM FROM THE ELBOW



With the EPITEC® FIX straps, pressure is applied through a viscoelastic cushion with special adjustment points. Correctly positioned, these reduce the force of traction which results from the insertion of the epicondyle musculature tendon; achieving this through the orientation of the muscular fibres.



ADJUSTABLE SILICONE

Silicone points which stick to the skin by pressure and impede slipping of the orthosis on the forearm.



FASTENING STRAP

It adjusts to the forearm and wraps around itself with a plastic buckle joined to an elastic regulating band which allows individualized adaptation of pressure according to a fixed scale. The area of closure is reinforced with soft Velcro to achieve a stronger and longer-lasting adjustment.



CUSHION

Made of viscoelastic material with special pressure points which apply a selective pressure at the same time as an effective massage and is lined on one side with a soft hypoallergenic towel and on the other with micro-velcro which allows it to stick to the orthosis and thus avoiding slipping.



SOUL OF THERMOPLASTIC

Molded and shaped to fit the shape of the upper arm and serves as subject to the locking pin buckle. When performing the compression positions of the forearm muscles causing decompression and relief.