

Star Brace® Dynamic Fix



Ref.: SD103

Ref.: SD104 HIGH



STAR BRACE® DYNAMIC FIX RIGID CORSET

Characteristics:

- Progressive lumbosacral orthosis, made of lightweight and breathable double-layer mesh material that facilitates ventilation and moisture removal.
- Pre-moulded thermoplastic lumbar module made of low-density polyethylene, which is secured with a textile cover that allows for removal and insertion.
- Lumbar plate that is interchangeable with the thermoplastic module for evolutionary immobilisation.
- Selective compression pulley system for the lumbar and sacral areas independently.
- The lacing pull tabs, with a thumbhole, facilitate donning and doffing the orthosis.
- Front Velcro strap that allows for the placement or removal of an abdominal plate to increase or decrease the level of containment if necessary.
- Includes an abdominal plate for greater containment in the area.

Indications:

- Disc hernias, osteoporosis, spondyloarthrosis, spondylolysis with or without stability, trauma, compression fractures, degenerative disc disease, postoperative fusion, chronic muscle disease, and treatment of vertebral metastases.

Colour: ● Black.

Available models:

- SD103: Short.
- SD104: High.

SIZES	PERIMETER CM	TROCANTER	FRONT HEIGHT	BACK HEIGHT
1	70-84			
2	84-98		16 cm (SD103) 18 cm (SD104)	25 cm (SD103) 30 cm (SD104)
3	98-112			
4	112-126			
5	126-140			

Components SD103/SD104




SUPPORT



THERMOPLASTIC LUMBAR MODULE



ABDOMINAL PLATE



LOW-DENSITY TP PLATE
SEMI-RIGID CONTAINMENT